

## Product Datasheet

### ATF4 (Phospho-Ser245) Antibody (orb1165412)

|                            |                                   |
|----------------------------|-----------------------------------|
| <b>Description</b>         | ATF4 (Phospho-Ser245)<br>Antibody |
| <b>Species/Host</b>        | Rabbit                            |
| <b>Reactivity</b>          | Human                             |
| <b>Conjugation</b>         | Unconjugated                      |
| <b>Tested Applications</b> | IHC, WB                           |
| <b>Note</b>                | For research use only             |
| <b>Clonality</b>           | Polyclonal                        |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |