

## Product Datasheet

### Phospho-Bad(S118) Antibody (orb1167079)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | Phospho-Bad(S118)<br>Antibody   |
| <b>Species/Host</b>        | Rabbit                          |
| <b>Reactivity</b>          | Human                           |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | IHC-P, WB                       |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Polyclonal                      |
| <b>Expiration Date</b>     | 12 months from date of receipt. |