

## Product Datasheet

### Phospho-ABL1(Y245) Antibody (orb1167087)

|                            |                                 |
|----------------------------|---------------------------------|
| <b>Description</b>         | Phospho-ABL1(Y245)<br>Antibody  |
| <b>Species/Host</b>        | Rabbit                          |
| <b>Reactivity</b>          | Human, Mouse                    |
| <b>Conjugation</b>         | Unconjugated                    |
| <b>Tested Applications</b> | IF, IHC-P, WB                   |
| <b>Note</b>                | For research use only           |
| <b>Clonality</b>           | Polyclonal                      |
| <b>Expiration Date</b>     | 12 months from date of receipt. |