


Product Datasheet

GAPDH antibody (orb1172451)

Description	GAPDH antibody
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant mouse GAPDH protein (1-333/333 aa)
Target	GAPDH
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P04406
Dilution Range	WB=1:5000-40000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human 293T cell; Lane...