

Product Datasheet

Calbindin D28K antibody (orb1172790)

Description	Calbindin D28K antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Target	Calbindin D28K
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1/Kappa
Clonality	Recombinant
Uniprot ID	P05937
Expiration Date	12 months from date of receipt.