

Product Datasheet

Acetyl-Histone H4 (Lys5/8/12/16) antibody (orb1172811)

Description	Acetyl-Histone H4 (Lys5/8/12/16) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Target	Acetyl-Histone H4 (Lys5/8/12/16)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
Uniprot ID	P62805
Expiration Date	12 months from date of receipt.