

Product Datasheet

Histone H3 (Phospho-Ser10) antibody (orb1172832)

Description	Histone H3 (Phospho-Ser10) antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Histone H3 around the phosphorylation site of Ser10 (RK(p-S)TG)
Target	Phospho-Histone H3 (Ser10)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1/Kappa
Clonality	Recombinant
MW	15 kDa
CAS Number	Histone H3.1
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.