

Product Datasheet

Mono-Methyl-Histone H3 (Lys79) antibody (orb1172860)

Description	Mono-Methyl-Histone H3 (Lys79) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Target	Mono-Methyl-Histone H3 (Lys79)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17 kDa
CAS Number	Histone H3.1
Uniprot ID	P68431
Expiration Date	12 months from date of receipt.