

Product Datasheet

COPS3 antibody (orb1173013)

Description	COPS3 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Target	COPS3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG/Kappa
Clonality	Recombinant
Uniprot ID	Q9UNS2
Expiration Date	12 months from date of receipt.