

Product Datasheet

COPE antibody (orb1173090)

Description	COPE antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, WB
Target	COPE
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1/Kappa
Clonality	Recombinant
MW	34 kDa
Uniprot ID	O14579
Expiration Date	12 months from date of receipt.