

Product Datasheet

CD5 antibody (orb1173094)

Description	CD5 antibody
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC
Target	CD5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2a/Kappa
Clonality	Recombinant
Uniprot ID	P06127
Expiration Date	12 months from date of receipt.