

Product Datasheet

Amino Acid (AA) Colorimetric Assay Kit (orb1173181)

Description	CheKine Micro Amino Acid (AA) Assay Kit is designed for detecting Amino Acid (AA) content in the sam...
Range	0.625-40 µmol/mL
Storage	Storage at 4°C and Keep from light immediately upon receipt. Kit has a storage time of 6 months from receipt. Refer to list of materials supplied for storage conditions of individual components.
Note	For research use only
Application notes	The kit provides a simple method for detecting FC concentration in a variety of biological samples such as serum, plasma, tissues, cells, and bacteria.
Sensitivity	0.5 µmol/mL
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	Please enquire.