

Product Datasheet

Glycogen Assay Kit (Colorimetric) (orb1173189)

Description	CheKine Micro Glycogen Assay Kit provides a simple, sensitive, rapid glycogen detection method, suit...
Storage	The kit is valid for 12 months. Please refer to the table below Materials supplied and Storage conditions to store all the components.
Note	For research use only
Application notes	<p>CheKine Micro Glycogen Assay kit provides a simple, convenient, rapid and colorimetric glycogen detection method. The principle of kit is strong alkaline extract glycogen, under strong acidic conditions, blue compounds with anthrone chromogenic agent with characteristic absorption peak at 620 nm. In a certain concentration range, glycogen concentration is linearly related to 620 nm absorbance. According to the standard curve, the glycogen concentration in the sample can be calculated.</p>
Hazard Information	<p>The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.</p>
Expiration Date	Please enquire.