

Product Datasheet

Glutathione S-Transferase (GST) Colorimetric Assay Kit (orb1173237)

Description	Biorbyt CheKine Micro Glutathione S-Transferase (GST) Assay Kit is designed for detecting GST activi...
Storage	Storage at 4°C and protected from light upon receipt. Kit has a storage time of 12 months form receipt.
Note	For research use only
Application notes	GST catalyzes the binding of GSH to CDNB. The conjugation is accompanied by an increase in absorbance at 340 nm. The rate of increase is directly proportional to the GST activity in the sample.
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	Please enquire.