

## Product Datasheet

### HRP Conjugated Anti-DDDDK Tag Mouse Monoclonal Antibody (1B10) (orb1173288)

## Description

## Featured HRP Conjugated Anti-

|                            |                                                                                                                                                                                                                                                                                                                                                                  |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                            | Antibody (1B10), specially designed for your...                                                                                                                                                                                                                                                                                                                  |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                                                                                                                                                                                            |
| <b>Conjugation</b>         | HRP                                                                                                                                                                                                                                                                                                                                                              |
| <b>Tested Applications</b> | WB                                                                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                                                                                                                                                                                                                |
| <b>Preservatives</b>       | Liquid in PBS, pH 7.4 and 50% Glycerol.                                                                                                                                                                                                                                                                                                                          |
| <b>Form/Appearance</b>     | Liquid solution                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>             | Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.                                                                                                                                                   |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                            |
| <b>Application notes</b>   | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:5000-1:10000).                                                                                                                                                                                                             |
| <b>Isotype</b>             | Mouse IgG                                                                                                                                                                                                                                                                                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                                                                                                                                                                                                       |
| <b>Hazard Information</b>  | The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                                                                                                  |