

Product Datasheet

Magnetic Beads Conjugated Anti-DDDDK Tag Mouse Monoclonal Antibody (1B10) (orb1173290)

Description
Featured Anti-DDDDK Tag Mouse

	Beads, specially designed for immunoprec...
Species/Host	Mouse
Conjugation	Magnetic Beads
Tested Applications	IP
Immunogen	Synthetic Peptide
Preservatives	Suspended in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative.
Form/Appearance	Liquid solution
Storage	Store at 4°C for 12 months. Avoid freeze-thaw or centrifugation.
Note	For research use only
Application notes	Optimal working dilutions should be determined experimentally by the investigator.
Isotype	Mouse IgG1
Clonality	Monoclonal
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road, Cambridge, CB5 8LA, United Kingdom

Email: info@biorbyt.com | Phone: +44 (0) 1223 859-353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

68 TW Alexander Drive
Research Triangle Park
Durham, North Carolina
27709, United States

Email: info@biorbyt.com | Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558