

Product Datasheet

GST-Tag Protein Purification Kit (Glutathione) (orb1173394)

Description	PurKine GST Tagged Protein Purification Protein Purification Kit
Storage	Stable for one year at 2-8°C from date of shipment. Do not freeze.
Note	For research use only
Application notes	PurKine GST-Tag Protein Purification Kit effectively purifies high levels of overexpressed GST-tagged fusion proteins in gravity column procedures at a variety of scales. The Resin consists of 90µm beads of cross- linked 4% agarose, to which reduced glutathione has been coupled. Tests confirm that performance equals or exceeds popular GST-Tag Purification Resins from other suppliers, and no decrease in performance occurs after at least five repeated uses.
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	Please enquire.