

Product Datasheet

APP-1/PABPC4 Antibody (orb1173479)

Description

APP-1/PABPC4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, ICC, IF, IHC, WB

Applications
Immunogen

A synthetic peptide corresponding to a sequence of human APP-1/PABPC4 (RPPYYTPNQLAQMRPNRWQQ).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW


75 kDa

Uniprot ID
Q13310
Dilution Range


Western blot, 0.25-0.5 µg/ml, Human, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Human Immunocytochemistry/Immunofluorescence, 5 µg/ml, Human Immunofluorescence, 5 µg/ml, Human Flow Cytometry, 1-3 µg/1x10⁶ cells, Human, Mouse

Expiration Date


12 months from date of receipt.



IHC analysis of APP-1/PABPC4 using anti-...



IHC analysis of APP-1/PABPC4 using anti-...



IF analysis of APP-1/PABPC4 using anti-A...