

Product Datasheet

Rat MBL1 ELISA Kit (orb1173622)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rat
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma and other biological fluids
Uniprot ID	P19999
Solubility (25°C)	10 ng/mL
Sensitivity	0.059 ng/mL
Expiration Date	Please enquire.

