

## Product Datasheet

### Zebrafish LEP ELISA Kit (orb1173623)

|                            |                                                                                                              |
|----------------------------|--------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...      |
| <b>Reactivity</b>          | Zebrafish                                                                                                    |
| <b>Range</b>               | 0.16-10 ng/mL                                                                                                |
| <b>Tested Applications</b> | ELISA                                                                                                        |
| <b>Note</b>                | For research use only                                                                                        |
| <b>Application notes</b>   | Sandwich                                                                                                     |
| <b>Sample Types</b>        | serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| <b>Solubility (25°C)</b>   | 10 ng/mL                                                                                                     |
| <b>Sensitivity</b>         | 0.058 ng/mL                                                                                                  |
| <b>Expiration Date</b>     | Please enquire.                                                                                              |

