

Product Datasheet

Mouse Chia ELISA Kit (orb1173624)

Description The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...

Reactivity Mouse

Range 0.16-10 ng/mL

Tested Applications ELISA

Note For research use only

Application notes Sandwich

Sample Types serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Uniprot ID [Q91XA9](#)

Solubility (25°C) 10 ng/mL

Sensitivity 0.058 ng/mL

Expiration Date Please enquire.

