


Product Datasheet


BAFF-R Antibody (orb1173654)

Description	Anti-BAFFR (ianalumab biosimilar) mAb
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	BAFF-R
Preservatives	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization.
Form/Appearance	Powder
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	IgG1
Clonality	Recombinant
Uniprot ID	Q96RJ3
Dilution Range	ELISA 1:5000-10000 Flow Cyt 1:100
Expiration Date	12 months from date of receipt.



ELISA plate
pre-coated by
2 µg/ml (100

...



HEK293 cell
line
transfected
with irrele...