

Product Datasheet

DNAJA1 Antibody (orb1174151)

Description	DNAJA1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, WB
Immunogen	A synthesized peptide derived from human DNAJA1
Target	DNAJA1
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Monoclonal
Uniprot ID	P31689
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200, FC: 1:50-1:200
Expiration Date	12 months from date of receipt.