



Product Datasheet

CD3 epsilon antibody (orb1175)


Description	CD3 epsilon antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse CD3 epsilon
Target	CD3E
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	P07766
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
Mouse
spleen...



Flow
cytometric
analysis of
(mo)
Thymocy...



Flow
cytometric
analysis of
(mo)
Thymocy...