

Product Datasheet

C1orf93 antibody (orb1175195)

Description	Rabbit polyclonal antibody to C1orf93
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 102-166aa of human C1orf93 protein
Target	PRXL2B
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	21kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.