

Product Datasheet

ACE2 antibody (orb1175405)

Description	Rabbit polyclonal antibody to ACE2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant fragment corresponding to 130-294aa of human ACE2 protein
Target	ACE2
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	805aa, 92KD; 555aa, 64KD
Dilution Range	ELISA: 1:20,000-1:80,000
Expiration Date	12 months from date of receipt.