

Product Datasheet

C1orf2 antibody (orb1175431)

Description	Rabbit polyclonal antibody to C1orf2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant fragment corresponding to 96-189aa of human ENTREP3 protein
Target	ENTREP3
Form/Appearance	Lyophilized powder with 1% rAlbumin and 0.02% NaN ₃ .
Storage	Shipped at ambient temperature or with ice bags if in liquid. Store immediately at -20°C upon receipt. Avoid repeating freeze-thaw cycles.
Note	For research use only
Application notes	Reconstitute in 100μL of sterile H ₂ O with 50% glycerol.
Isotype	IgG
Clonality	Polyclonal
MW	69kDa
Dilution Range	WB: 1:200-1:1,000
Expiration Date	12 months from date of receipt.