

Product Datasheet

Monoclonal Anti-IFN γ antibody, Human (8C5F8) (orb1176541)

Description	Monoclonal Anti-IFN γ antibody, Human (8C5F8)
Reactivity	Human
Endotoxins	1.0 EU per μ g
Conjugation	Unconjugated
Target	IFN-gamma
Preservatives	Powder
Form/Appearance	Powder
Storage	-20 to -70°C for 12 months in lyophilized state from date of receipt. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Note	For research use only
Application notes	This antibody can be paired with other IFN γ antibodies to detect IFN γ in sandwich ELISA or LFA assay.
Isotype	Human IgG1 Human Kappa
Clonality	Recombinant
Purity	95%
Expiration Date	12 months from date of receipt.