

Product Datasheet

Human LTBR / TNFRSF3 Protein (orb1176544)

Description

Human LTBR, Fc Tag is expressed from

Gln 31 - Met 227 (Acce...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

LTBR

Preservatives

50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 48.8 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_002333.1

Purity

95%

Source

Human LTBR, Fc Tag (orb1176544) is expressed from human 293 cells (HEK293). It contains AA Gln 31 - Met 227 (Accession # P36941-1).

MW

48.8 kDa

Uniprot ID
P36941
NCBI
NP_002333.1
Expiration Date

6 months from date of receipt.