

Product Datasheet

Human IL-3 R Beta / CD131 Protein (orb1176553)

Description CD131- the common β Subunit, is a member of the type 1 cytokine receptor family. There are three ma...

Reactivity Human

Endotoxins 1.0 EU per μ g

Conjugation Unconjugated

Target CD131

Preservatives 51 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

Form/Appearance Powder

Storage -20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 75.1 kDa. The protein migrates as 80-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.


Purity 95%

Source Human CD131, Fc Tag (orb1176553) is expressed from human 293 cells (HEK293). It contains AA Trp 17 - Trp 443 (Accession # P32927-1).

MW 75.1 kDa

NCBI [NP_000386.1](#)

Expiration Date 6 months from date of receipt.



Loaded
Human
CD131, Fc
Tag (Cat.
No. orb...