

Product Datasheet

Mouse Catenin Beta CTNNb1 ELISA Kit (orb1176636)

Description	Mouse Catenin Beta CTNNb1 ELISA Kit
Range	250-20,000 pg/mL
Tested Applications	ELISA
Target	Mouse Catenin Beta CTNNb1
Storage	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
Note	For research use only
Sample Types	100 µL final volume, dilution factor varies on samples
Sensitivity	100 pg/mL
Expiration Date	Please enquire.