

Product Datasheet

Mouse Betacellulin BTC ELISA Kit (orb1176639)

Description	Mouse Betacellulin BTC ELISA Kit
Range	7.8 - 500 pg/mL
Tested Applications	ELISA
Target	Mouse Betacellulin BTC
Storage	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
Note	For research use only
Sample Types	100 µL final volume, dilution factor varies on samples
Sensitivity	10.0 pg/mL
Expiration Date	Please enquire.