

## Product Datasheet

### Dog/Canine Interleukin 8 IL-8 ELISA Kit (orb1176672)

|                            |                                                                          |
|----------------------------|--------------------------------------------------------------------------|
| <b>Description</b>         | Dog/Canine Interleukin 8 IL-8 ELISA Kit                                  |
| <b>Range</b>               | 15.6-1000 pg/mL                                                          |
| <b>Tested Applications</b> | ELISA                                                                    |
| <b>Target</b>              | Dog/Canine Interleukin 8 IL-8                                            |
| <b>Storage</b>             | Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C |
| <b>Note</b>                | For research use only                                                    |
| <b>Sample Types</b>        | 100 µL final volume, dilution factor varies on samples                   |
| <b>Sensitivity</b>         | 1.0 pg/mL                                                                |
| <b>Expiration Date</b>     | Please enquire.                                                          |