

Product Datasheet

Mouse Fibrinogen (FBG) ELISA Kit (orb1176734)

Description	Mouse Fibrinogen (FBG) ELISA Kit
Tested Applications	ELISA
Target	Mouse Fibrinogen (FBG)
Storage	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
Note	For research use only
Sample Types	50 µL final volume, dilution factor varies on samples
Sensitivity	3.0 ng/ml
Expiration Date	Please enquire.