

Product Datasheet

Human CD19 SAFIRE antibody (orb1176810)

Description	Mouse monoclonal antibody to CD19 SAFIRE
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Target	Human CD19 SAFIRE
Preservatives	PBS (pH 7.2) with 0.09% sodium azide. May contain a stabilizer or a carrier protein
Form/Appearance	Lyophilized
Storage	Store at 4 °C, as supplied. Do not freeze.
Note	For research use only
Isotype	IgG1, kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.