

Product Datasheet

Cathepsin D antibody (orb1176828)

Description	Mouse monoclonal antibody to CTSD
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Human liver cathepsin D.
Target	CTSD
Preservatives	Mouse ascites fluid, 1.2% sodium acetate, 2 mg rAlbumin, with 0.01 mg NaN ₃ as preservative. *carrier free antibody available upon request
Form/Appearance	Lyophilized
Storage	At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	Western blot: 0.5-1µg/ml. Immunohistochemistry(Paraffin-embedded Section): 1-2µg/ml
Expiration Date	12 months from date of receipt.