

## Product Datasheet

### beta-Amyloid Protein antibody (orb1176831)

|                            |                                                                                                                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse monoclonal antibody to beta-Amyloid Protein                                                                                                                                    |
| <b>Species/Host</b>        | Mouse                                                                                                                                                                                |
| <b>Reactivity</b>          | Human                                                                                                                                                                                |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                         |
| <b>Tested Applications</b> | IHC-P                                                                                                                                                                                |
| <b>Target</b>              | beta-Amyloid Protein                                                                                                                                                                 |
| <b>Preservatives</b>       | Mouse ascites fluid, 1.2% sodium acetate, 2 mg rAlbumin, with 0.01 mg NaN <sub>3</sub> as preservative. *carrier free antibody available upon request.                               |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                          |
| <b>Storage</b>             | At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing. |
| <b>Note</b>                | For research use only                                                                                                                                                                |
| <b>Isotype</b>             | IgG                                                                                                                                                                                  |
| <b>Clonality</b>           | Monoclonal                                                                                                                                                                           |
| <b>Dilution Range</b>      | Immunohistochemistry(Paraffin-embedded Section): 0.5-1µg/ml                                                                                                                          |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                      |