

Product Datasheet

5HT2A Receptor antibody (orb1177025)

Description	Rabbit polyclonal antibody to HTR2A
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human 5HT2A Receptor (400-431aa KENKKPLQLILVNTIPALAYKSSQLQMGQKKN) , different from the related mouse sequence by three amino acids, and from the related rat sequence by two amino acids.
Target	HTR2A
Preservatives	Antibody is lyophilized with 5 mg rAlbumin, 0.9 mg NaCl, 0.2 mg Na2HPO4, 0.05 mg NaN3. *carrier free antibody available upon request
Form/Appearance	Lyophilized
Storage	At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Western blot: 0.1-0.5µg/ml
Expiration Date	12 months from date of receipt.