

Product Datasheet

SFTPA1 antibody (orb1177358)

Description	Rabbit polyclonal antibody to SFTPA1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	SFTPA1
Preservatives	Antibody is lyophilized with 50% glycerol, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ . *carrier free antibody available upon request.
Form/Appearance	Liquid
Storage	At -20 °C for 12 months, as supplied, at 2-8 °C for one month. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Western blot: 0.1-0.5µg/ml
Expiration Date	12 months from date of receipt.