

Product Datasheet

gamma-Tubulin antibody (orb1178262)

Description	Mouse monoclonal antibody to GGA1; GGA2; GGA3
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Target	GGA1; GGA2; GGA3
Preservatives	Mouse ascites fluid, 1.2% sodium acetate, 2 mg rAlbumin, with 0.01 mg NaN ₃ as preservative. *carrier free antibody available upon request.
Form/Appearance	Lyophilized
Storage	At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Dilution Range	Western blot: 0.25-0.5µg/ml. Immunocytochemistry: Suitable
Expiration Date	12 months from date of receipt.