

Product Datasheet

Dopamine Receptor D1 antibody (orb1178421)

Description	Rabbit polyclonal antibody to DRD1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	DRD1
Preservatives	Antibody is lyophilized with 5 mg rAlbumin, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ , 0.05 mg Thimerosal and 0.05 mg NaN ₃ . *carrier free antibody available upon request.
Form/Appearance	Lyophilized
Storage	At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Western blot: 0.1-0.5µg/ml
Expiration Date	12 months from date of receipt.