

## Product Datasheet

### Chk2 antibody (orb1178571)

|                            |                                                                                                                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CHEK2                                                                                                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                               |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                                    |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                         |
| <b>Tested Applications</b> | WB                                                                                                                                                                                   |
| <b>Target</b>              | CHEK2                                                                                                                                                                                |
| <b>Preservatives</b>       | Antibody is lyophilized with 5 mg rAlbumin, 0.9 mg NaCl, 0.2 mg Na2HPO4, 0.05 mg Thimerosal and 0.05 mg NaN3. *carrier free antibody available upon request.                         |
| <b>Form/Appearance</b>     | Lyophilized                                                                                                                                                                          |
| <b>Storage</b>             | At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing. |
| <b>Note</b>                | For research use only                                                                                                                                                                |
| <b>Isotype</b>             | IgG                                                                                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                           |
| <b>Dilution Range</b>      | Western blot: 0.1-0.5µg/ml                                                                                                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                      |