

Product Datasheet

Mouse anti-cMyc tag mAb (orb1178846)

Description	Mouse anti-cMyc tag mAb is specific to EQKLISEEDL tags which can be applied to several immunoassays,...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IF, IP, WB
Immunogen	Peptide: EQKLISEEDL (cMyc)
Preservatives	Phosphate Buffered Saline containing 0.01% thimerosal, pH 7.4.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Monoclonal
Purity	> 95% (SDS-PAGE)
Expiration Date	12 months from date of receipt.