

Product Datasheet

SARS-CoV-2 Nucleocapsid Protein (Omicron B1.1.529 Variant), Tag Free, HEK293 (orb1178856)

Description	SARS-CoV-2 Nucleocapsid Protein (Omicron B1.1.529 Variant), Tag Free, HEK293
Form/Appearance	Lyophilized
Storage	Lyophilized protein should be stored at -20°C. Proteins are stable for up to one year from date of receipt at -20°C. Upon reconstitution, protein aliquots should be stored at -20°C or -80°C. Please avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	It is recommended to reconstitute the lyophilized protein in sterile H ₂ O to a concentration not less than 100 µg/mL and incubate the stock solution for at least 20 min to ensure sufficient re-dissolved.
Expiration Date	6 months from date of receipt.