

Product Datasheet

SARS-CoV-2 Trimeric Spike Protein (Beta B1.351 Variant), His Tag, HEK293 (orb1178867)

Description	SARS-CoV-2 Trimeric Spike Protein (Beta B1.351 Variant), His Tag, HEK293
Form/Appearance	Lyophilized
Storage	Stored at -80°C. This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles.
Note	For research use only
Application notes	Liquid, Phosphate buffered saline, pH 7.4.
Expiration Date	6 months from date of receipt.