

Product Datasheet

Human anti-SARS-CoV & CoV-2 Spike antibody (IgM), clone CR3022 (orb1178881)

Description	Human anti-SARS-CoV-2 Spike Antibody [CR3022] recognizes human SARS-CoV and CoV-2 Spike protein with...
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	Phosphate Buffered Saline, pH 7.4.
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store at 4°C for six months.
Note	For research use only
Isotype	IgM
Clonality	Monoclonal
Dilution Range	ELISA 1:5000-20000
Expiration Date	12 months from date of receipt.