

Product Datasheet

A 784168 (orb1179313)

Description	A 784168
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{19}H_{15}F_6N_3O_3S$
Purity	> 98%
MW	479.4
CAS Number	[824982-41-4]
SMILES	<chem>FC(F)(F)C1=CC=CN=C1N2CCC(=CC2)C(=O)NC3=CC=C(C=C3)S(=O)(=O)C(F)(F)F</chem>
Expiration Date	12 months from date of receipt.