

Product Datasheet

SAFit1 (orb1180490)

Description	SAFit1
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₄₂ H ₅₃ NO ₁₁
Purity	> 98%
MW	747.87
CAS Number	[1643125-32-9]
SMILES	<chem>O=C(N1[C@@H](CCCC1)C(O[C@@H](C2=CC(OCC(O)=O)=CC=C2)CCC3=CC(OC)=C(OC)C=C3)=O)[C@@H](C4CCCC4)C5=CC(OC)=C(OC)C(OC)=C5</chem>
Expiration Date	12 months from date of receipt.