

Product Datasheet

QCA570 (orb1180496)

Description	QCA570
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	$C_{39}H_{33}N_7O_4S$
Purity	> 98%
MW	695.79
CAS Number	[2207569-08-0]
SMILES	<chem>O=C(N(C1CCC(NC1=O)=O)C2)C3=C2C(C#CCCCN4N=CC(C#CC(S5)=C(CC6=CC=CC=C6)C(COC7)=C5N8C7=NN=C8C)=C4)=CC=C3</chem>
Expiration Date	12 months from date of receipt.