

Product Datasheet

YAP/TAZ inhibitor-1 (orb1180498)

Description	YAP/TAZ inhibitor-1
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₃₃ H ₃₉ N ₃ O ₅ S ₂
Purity	> 98%
MW	621.81
CAS Number	[2093565-23-0]
SMILES	<chem>O=S(C1=CC2=C(C=C1)C3=C(C=C(S(=O)(C4CCC(C)(C)CC4)=O)C=C3)/C2=NO)(N5CCN(C(C6=CC=CC=C6)C)CC5)=O</chem>
Expiration Date	12 months from date of receipt.